



Lobby dust-pan - Yellow

Lobby dust-pan to safely collect even large quantities of waste.



Code	Colour	Size (mm)	Material	Net weight (kg)
1055Y	Yellow	275 x 310 x 950	PP - AL	0,783

Usage information

Minimum and maximum temperature	min -20 / max 100 C°
Washing (dishwasher)	-
Autoclave	121° C for 20 min
User and maintenance manual	- Remove packaging and any label from tools and disinfect them. - Both manual and autoclave disinfection is possible. - Sanitize tools using appropriate chemical decontamination products, in the concentrations, times and temperatures provided. - Make sure tools completely dry up before next use. - Hang tools while drying up to avoid any contamination or breakage. - The storing area must be clean and far away from sources of intense light or heat. - Do not use damaged tools.

Compliance

- FDA

Logistical information

Internal packaging type	-	Product barcode	8008805155507
External packaging type	Carton	Box barcode	08008805155507
Quantity per pack (pcs)	6	Customs tariff	96039099
Box gross weight (kg)	5,496	Country of origin	Mexico
Box dimensions (mm)	870 x 320 x 320	Palletising (pcs)	144
Box volume (m³)	0,089088	Pallet dimensions (cm)	130 x 90 x 206,5

Available versions

Size		Material	
Colour	B R G W		